

Gyro™

Fixed speed microcentrifuge
QP-1800-01

TECHNICAL SPECIFICATIONS



FEATURES

Single speed microcentrifuge, robust and user-friendly, with interchangeable rotors.

SPEED

10,000 RPM

MAXIMUM RCF

4,800 g

ROTOR CAPACITY

6 x 1.5 ml, 16 x 0.2 ml, 6 adapters for 0.2 or 0.5 ml

NOISE LEVEL

< 52 dB

DIMENSIONS

17 x 15 x 13 cm

WEIGHT

1 kg

VOLTAGE

AC 100-240 V, 50/60 Hz

ADDITIONAL FEATURES

Automatic safety stop

1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com