

P51 Pro™ Molecular Fluorescence Viewer



QP-1900-10

TECHNICAL SPECIFICATIONS

ILLUMINATION	HIGH INTENSITY BLUE LEDs, ~465 NM PEAK WAVELENGTH
FILTER	AMBER FOR BLUE LIGHT FILTERING
INTENSITY	VARIABLE
POWER	5 V @ 0.5 A USB-C OR PD POWER PACK
RACK CAPACITY	8 POSITIONS
TUBE FORMAT	0.2 ML TUBES, DOMED OR FLAT CAPS, STRIP COMPATIBLE
DIMENSIONS	4.1" x 2.5" x 6.7" (10.4 x 6.3 x 17 CM)
WEIGHT	0.56 LB (255 G)

1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com